

G-Force™ Series

STRETCHER SHEETS




**Color Coded
Lift Capacity™**

Taylor's Color Coded Lift Capacity™ fitted stretcher sheets set the new industry standard for transferring a patient with confidence. By utilizing Taylor's Color Coded Lift Capacity™ system you may easily choose the Taylor stretcher sheets best suited for your EMS service. The Taylor G-Force™ Fitted Sheet Series is designed and tested in the transfer of patients up to 450 lbs. The Taylor G-Force™ Fitted Stretcher Sheet Series has you covered for both lift capacity and price point.

- The Color Coded sheets are designed to aid in patient transfers for the following weights when using proper lifting method:
 - Green – 250 lbs.
 - Blue – 400 lbs.
 - White – 450 lbs.

- Custom fitted to the unique Stryker® XPS and Ferno® Power X1 mattress designs
- Available in fluid resistant*, fluid impervious** and heavy duty versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets



Item#	Product	Color	Material	Size	Packaging
95-ECOFORCE	ECO-FORCE™ Fitted Sheet	Green	Fluid Resistant	36" x 90"	50/case
92-EFORCE	EFORCE™ Fitted Sheet	Grey	Fluid Resistant	36" x 90"	50/case
90-GFRC3690	G-Force™ Fitted Sheet	Blue	Fluid Resistant	36" x 90"	30/case
90-IFRC3690	I-Force™ Fitted Sheet	Blue	Fluid Impervious	36" x 90"	30/case
90-SUPER-G	SUPER-G™ Fitted Sheet	White	Fluid Resistant	36" x 90"	30/case

*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure
 **Fluid Impervious - Fluid does not penetrate the material even with time and pressure

