

# DISPOSABLE FITTED STRETCHER SHEETS

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's durable fitted stretcher sheets are the best disposable fitted sheets available! Specially designed for patient transfer, Taylor fitted sheets provide the surest fit, gripping the mattress tightly and preventing slippage.

- The Color Coded sheets are designed to aid in patient transfers for the following weights when using proper lifting method:

Green – 250 lbs.      Blue – 400 lbs.  
Grey – 300 lbs.      White – 450 lbs.


- Available in fluid resistant\* and fluid impervious\*\* versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets
- Our selection includes fitted sheets designed to aid in patient transfers up to 450 lbs.
- Options available for standard and large cots, stretchers & gurneys



Taylor™  
G-Force Series




Taylor™  
SureFit Series



Item#	Product	Color	Material	Size	Packaging
<b>G-FORCE SERIES</b>					
95-ECOFORCE	ECOFORCE™ Fitted Sheet for G-FORCE™ shaped mattresses	Green	Fluid Resistant	36" x 90"	50/case
92-EFORCE	EFORCE™ Fitted Sheet for G-FORCE™ shaped mattresses	Grey	Fluid Resistant	36" x 90"	50/case
90-GFRC3690	G-Force™ Fitted Sheet for G-FORCE™ shaped mattresses	Blue	Fluid Resistant	36" x 90"	30/case
90-IFRC3690	I-Force™ Fitted Sheet for G-FORCE™ shaped mattresses	Blue	Fluid Impervious	36" x 90"	30/case
90-SUPER-G	Super-G™ Fitted Sheet for G-FORCE™ shaped mattresses	White	Fluid Resistant	36" x 90"	30/case
<b>SUREFIT SERIES</b>					
90-BSS3084	Standard Fitted Sheet 30"	Blue	Fluid Resistant	30" x 84"	50/case
90-BIS3084	Impervious Fitted Sheet 30"	Blue	Fluid Impervious	30" x 84"	50/case
90-BSS3484	Standard Fitted Sheet 34"	Blue	Fluid Resistant	34" x 84"	50/case
90-BIS3484	Impervious Fitted Sheet 34"	Blue	Fluid Impervious	34" x 84"	50/case
90-BXS3490	Extended Length Fitted Sheet	Blue	Fluid Resistant	34" x 90"	50/case
90-IXS3490	Impervious Extended Length Fitted Sheet	Blue	Fluid Impervious	34" x 90"	50/case
90-WHS3484	Heavy Duty Fitted Sheet	White	Fluid Resistant	34" x 84"	50/case
92-SFS3084	Standard Fitted Sheet	Grey	Fluid Resistant	30" x 84"	50/case
92-IFS3084	Impervious Fitted Sheet	Grey	Fluid Impervious	30" x 84"	50/case
95-EFS3084	Economy Fitted Sheet	Green	Fluid Resistant	30" x 84"	50/case

\*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure

\*\*Fluid Impervious - Fluid does not penetrate the material even with time and pressure

