


TITAN™ & TITAN II™

SOFT STRETCHERS

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service! The Titan Series are products no service can do without. They are essential to keep patients and rescuers safe as they navigate challenges faced daily.

- Cost effective, reusable, easy to store and durable
- The essential tool to aid in lifting bariatric patients
- Vital to have on the scene when space limitations prohibit use of standard EMS stretchers
- 16 SureGrip™ rubber handles ensure rescuers the comfort and stability of a firm grip, which provides easy lifting while transporting a patient
- Made of ultra-strong PVC and nylon webbing
- TITAN is designed to accommodate 16" – 18" backboards*
- TITAN II is designed to be used without a backboard
- Easily disinfected and air dried
- Folds for space-saving storage in the accompanying custom Titan carrying case
- Available individually or in cases of 5
- Lifting capacity of 1500 lbs



| Item# | Product | Color | Material | Size | Packaging |
|-------------|----------|-------|--|-----------|------------|
| 44-TT4080 | Titan | Black | PVC, nylon webbing, SureGrip™ rubber handles | 40" x 80" | 5 per case |
| 44-TITAN II | Titan II | Black | PVC, nylon webbing, SureGrip™ rubber handles | 40" x 80" | 5 per case |

*Backboard not included



LP44A0823