


TITAN™

SOFT STRETCHERS

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service! The Taylor Titan is the product no service can do without. The Titan Soft Stretcher is essential to keep patients and rescuers safe as they navigate challenges faced daily.

- Cost effective, reusable, easy to store and durable
- The essential tool to aid in lifting bariatric patients
- Vital to have on the scene when space limitations prohibit use of standard EMS stretchers
- 16 SureGrip™ rubber handles ensure rescuers the comfort and stability of a firm grip, which provides easy lifting while transporting a patient
- Made of ultra-strong PVC and nylon webbing
- Designed for use with or without a backboard. Accommodates 16" – 18" backboards
- Easily disinfected and air dried
- Folds for space-saving storage in the accompanying custom Titan carrying case
- Available individually or in cases of 5
- Lifting capacity of 1500 lbs



Item #	Product	Color	Material	Size	Packaging
44-TT4080	Titan™	Black	PVC, nylon webbing, SureGrip™ rubber handles	40" x 80"	5 per case

*Backboard not included



LP44A0823