



G-Force™ Series

STRETCHER SHEETS




Color Coded
Lift Capacity™

Taylor's Color Coded Lift Capacity™ fitted stretcher sheets set the new industry standard for transferring a patient with confidence. By utilizing Taylor's Color Coded Lift Capacity™ system you may easily choose the Taylor stretcher sheets best suited for your EMS service. The Taylor G-Force™ Fitted Sheet Series is designed and tested in the transfer of patients up to 450 lbs. The Taylor G-Force™ Fitted Stretcher Sheet Series has you covered for both lift capacity and price point.

- The Color Coded sheets are designed to aid in patient transfers for the following weights when using proper lifting method:
 - Green – 250 lbs.
 - Blue – 400 lbs.
 - White – 450 lbs.

- Custom fitted to the unique Stryker® XPS and Ferno® Power X1 mattress designs
- Available in fluid resistant*, fluid impervious** and heavy duty versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets



| Item# | Product | Color | Material | Size | Packaging |
|-------------|-------------------------|-------|------------------|-----------|-----------|
| 95-ECOFORCE | ECO-FORCE™ Fitted Sheet | Green | Fluid Resistant | 36" x 90" | 50/case |
| 92-EFORCE | EFORCE™ Fitted Sheet | Grey | Fluid Resistant | 36" x 90" | 50/case |
| 90-GFRC3690 | G-Force™ Fitted Sheet | Blue | Fluid Resistant | 36" x 90" | 30/case |
| 90-IFRC3690 | I-Force™ Fitted Sheet | Blue | Fluid Impervious | 36" x 90" | 30/case |
| 90-SUPER-G | SUPER-G™ Fitted Sheet | White | Fluid Resistant | 36" x 90" | 30/case |

*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure
 **Fluid Impervious - Fluid does not penetrate the material even with time and pressure



LP90GFCC0823