

Mattress Covers


Taylor Healthcare Products, Inc. continues to be the industry leader in quality, innovation and customer service!

Through unsparing conditions and relentless transports, Taylor Healthcare has you covered.

Our durable, cost effective Mattress Cover protects and aids in extending the life of your cot mattress.

The SureFit™ Mattress Cover features a pocketed zipper design on one end, while the G-Force™ Mattress Cover utilizes an underside zipper placement.

- Availability for both SureFit™ and G-Force™ shape mattresses
- Easily cleaned and disinfected – featuring a smooth surface that saves time from having to clean in the joints and folds of the mattress
- Fluid impervious* material provides an extra barrier of protection
- Reduces risk of cross-contamination and improves infection control protocols



| Item# | Product | Color | Material | Size | Packaging |
|------------|-------------------------|-------|----------|-----------|------------|
| 46-SUREFIT | SureFit™ Mattress Cover | Black | PVC | 22" x 72" | 5 per Case |
| 46-GFORCE | G-FORCE™ Mattress Cover | Black | PVC | 26" x 74" | 5 per Case |

*Fluid Impervious - Fluid does not penetrate the material even with time and pressure



LP46MC1223