



G-Force™ Mattress Cover



Taylor Healthcare Products, Inc.

continues to be the industry
leader in quality, innovation
and customer service!

Through unsparing conditions
and relentless transports,
Taylor Healthcare has you covered.

Our durable, cost effective

Mattress Cover protects and aids in
extending the life of your cot mattress.

The G-Force™ Mattress Cover has an
underside zipper placement.

- Easily cleaned and disinfected, featuring a smooth surface that saves time from having to clean in to the joints and folds of the mattress
- Fluid impervious* material provides an extra barrier of protection
- Reduces risk of cross-contamination and improves infection control protocols



Item #	Product	Color	Material	Size	Packaging
46-GFORCE	Taylor G-Force™ Mattress Cover for G-Force shaped cots	Black	PVC	26" X 74"	5/case
46-GFORCE-G	Taylor G-Force™ Mattress Cover for G-Force shaped cots	Grey	PVC	26" X 74"	5/case

*Fluid Impervious- Fluid does not penetrate the material even with time and pressure with the exception of the zipper