

TITAN

Rapid Cooling Bag

Taylor Healthcare Products, Inc. continues to be the industry leader in quality, innovation and customer service!

Taylor Titan Rapid Cooling Bags set the new standard with efficiency, containment, functionality and strength for lifesaving solutions in heat related situations. Our bags consist of a six-handle design for stability while lifting and transporting patients.

- The cooling bags are designed and tested to aid in patient transport up to 575 lbs.

- Fluid Impervious*
- Made of ultra-strong PVC for durability



Front



Back

Item #	Product	Color	Material	Size	Packaging
46-TAY400460	Titan™ Rapid Cooling Bag	Royal Blue	PVC	44" x 95"	8/case

* Fluid Impervious - Fluid does not penetrate the material even with time and pressure with the exception of the zipper

**Photos display examples of use



LP46BBJ0624