

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's Storage Trays attach to stretchers and provide a convenient, easily accessible area to place medical bags, supplies or other small, often used equipment. Storage trays keep your hands free to maneuver the stretcher, maintain your back's health and improve patient care.

- Adds convenient storage space
- Constructed of durable nylon
- Secures easily to the stretcher with strong Velcro® fasteners
- Straps help secure equipment
- Head end trays come flat or with four storage pouches
- Fits most Stryker® cots




Headrest 4-Pouch Tray



Headrest 2-Pouch Tray



Head End Flat Tray



Headrest 3-Pouch Tray



Lower Bottom Flat Tray

Item#	Product	Color	Material	Size	Packaging
45-BHSP1035	Headrest 2-Pouch Tray	Black	Nylon	10" x 12"	10 per case
45-BHSP1030	Headrest 3-Pouch Tray	Black	Nylon	10" x 12"	10 per case
45-BHT1010P	Headrest 4-Pouch Tray	Black	Nylon	21" x 20"	10 per case
45-BHT1010F	Head End Flat Tray	Black	Nylon	21" x 20"	10 per case
45-BHT1020	Head End Retracting Flat Tray	Black	Nylon	21" x 20"	10 per case
45-BBT2010	Lower Bottom Flat Tray	Black	Nylon	38" x 18"	10 per case

*Misc. items displayed in the pouches are not included

LP45ST030425