

I-Force™ 3-Layer Fitted Stretcher Sheet

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Introducing our I-Force™ 3-Layer Fitted Stretcher Sheet – a specialty design meeting a need in the EMS industry. Reduce cross-contamination and debris transfer with the added feature of the impervious* foot pouch.

The three layers include a super absorbent top layer which quickly absorbs fluids, nonwoven layer for strength and the impervious poly layer prevents fluid from penetrating the material even with time and pressure.



- Aid in patient transfers up to 400 lbs. when using the proper lifting method

- Fluid Impervious* sheet and foot pouch

- Helps protect the mattress and patient against cross-contamination

- Custom fitted to the unique mattress designs of Ferno® Power X1 and X2; Stryker® XPS and Power-PRO

- Impervious foot pouch reduces debris transfer

- Elastic fitted design eliminates sheet slippage

- Soft feel maximizes patient comfort

- Helps keep responders' uniforms clean



Item #	Product	Color	Details	Size	Packaging
90-IFRC3690-3L	I-Force™ 3-Layer Fitted Stretcher Sheet	Blue/White	Fluid Impervious	36" x 90"	20 per case

*Fluid Impervious - Fluid does not penetrate the material even with time and pressure

