

The new standard
in
Cold Water
Immersion

RAPID COOLING BAG



Item#

48-TAY460KT

Features of the 48-TAY460KT

- Includes a Rapid Cooling Bag with 10 handles, 1 access port on front and 2 drain ports on back, 1 drain hose, double welded seams and 2 zipper pulls
- The Rapid Cooling Bag has a 575 lb. lift capacity
- Includes convenient Carry Bag
- Includes TitanTS™ soft stretcher with 12 handles and a 500 lb. lift capacity (Part# 43-TITANTS)

Item#

48-TAY460KT-PC

Features of the 48-TAY460KT-PC


- Includes a Rapid Cooling Bag with 10 handles, 1 access port on front and 2 drain ports on back, 1 drain hose, double welded seams and 2 zipper pulls
- The Rapid Cooling Bag has a 575 lb. lift capacity
- Includes convenient Carry Bag
- Includes TitanTS™ soft stretcher with 12 handles and a 500 lb. lift capacity (Part# 43-TITANTS)
- Includes a TitanPC™ soft stretcher with 14 handles and a 1,500 lb. lift capacity (Part# 43-TITANPC)



Item #	Product	Color	Material	Size	Packaging
46-TAY46ORCB-26	Titan™ Rapid Cooling Bag includes TitanTS™ Soft Stretcher (Black)	Royal Blue	PVC	44" x 95"	5/case
40-TAY46OCB	Carry Bag for Titan™ Rapid Cooling Bag	Black	PVC	22" x 5" x 15"	6/case
48-TAY460KT	Titan™ Rapid Cooling Bag Kit includes TitanTS™ & Carry Bag	Royal Blue	PVC	44" x 95"	Individual
48-TAY460KT-PC	Titan™ Rapid Cooling Bag Kit includes TitanTS™, TitanPC™ & Carry Bag	Royal Blue	PVC	44" x 95"	Individual
43-TITANPC	TitanPC™ Soft Stretcher	White	Nonwoven/Poly	40" x 80"	10/case
43-TITANTS	TitanTS™ Soft Stretcher	Black	Nonwoven/Poly	30" x 72"	20/case

The new standard
in
Cold Water
Immersion

RAPID COOLING BAG



Front View


Back View




Handles



Drain Hose



Front Access Port



Drain Ports on back

Item#

46-TAY46ORCB-26

Features of the 46-TAY46ORCB-26

- 10 handles for stability
- 1 access port on front and 2 drain ports on back
- Convenient drain hose
- Double welded seams ensure fluid containment
- Dual easy grip zipper pulls
- Lift capacity 575 lbs.
- Includes TitanTS™ soft stretcher with 12 handles and a 500 lb. lift capacity (Part# 43-TITANTS)

More to see



Item #	Product	Color	Material	Size	Packaging
46-TAY46ORCB-26	Titan™ Rapid Cooling Bag includes TitanTS™ Soft Stretcher (Black)	Royal Blue	PVC	44" x 95"	5/case
40-TAY46OCB	Carry Bag for Titan™ Rapid Cooling Bag	Black	PVC	22" x 5" x 15"	6/case
48-TAY46OKT	Titan™ Rapid Cooling Bag Kit includes TitanTS™ & Carry Bag	Royal Blue	PVC	44" x 95"	Individual
48-TAY46OKT-PC	Titan™ Rapid Cooling Bag Kit includes TitanTS™, TitanPC™ & Carry Bag	Royal Blue	PVC	44" x 95"	Individual
43-TITANPC	TitanPC™ Soft Stretcher	White	Nonwoven/Poly	40" x 80"	10/case
43-TITANTS	TitanTS™ Soft Stretcher	Black	Nonwoven/Poly	30" x 72"	20/case