


TITAN^{II} Soft Stretcher

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

The Taylor TitanII™ is the product no service can do without. The TitanII™ Soft Stretcher is essential to keeping patients and rescuers safe as they navigate challenges faced daily.

- Cost effective, reusable, easy to store and durable
- The essential tool to aid in lifting bariatric patients
- Vital to have on the scene when space limitations prohibit use of standard EMS stretchers
- 16 SureGrip™ rubber handles ensure rescuers the comfort and stability of a firm grip, which provides easy lifting while transporting a patient
- New enhancement features rivets throughout to secure handles
- Made of ultra-strong PVC and nylon webbing
- Easily disinfected and air dried
- Folds for space-saving storage in the accompanying custom TitanII™ carrying case
- Available individually or in cases of 5
- Lifting capacity of 1500 lbs.



Item #	Product	Color	Material	Size	Packaging
44-TITANII	TitanII™ Soft Stretcher	Black	PVC, nylon webbing, SureGrip™ rubber handles	40" x 80"	5/case

