

TITAN™ PC SERIES

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Our new TitanPC™, TitanPCX™ and TitanChair™ are beneficial additions to our Titan™ line of soft stretchers to aid in patient rescues and transfers when EMS professionals need compact and lightweight transport units to traverse confined spaces.

- Limited Use Soft Stretchers
- Strong, compact, lightweight and cost effective
- Each unit is easy to store, durable and fluid impervious*
- Vital to have on the scene when space limitations prohibit use of standard EMS stretchers
- **TitanPC™** and **TitanPCX™** each have 14 handles to ensure easy lifting while transporting a patient. Both designed to securely lift 1,500 lbs.
- The **TitanPCX™** accommodates 16"-18" backboards**
- The **TitanChair™** has 8 handles allowing responders the agility to transport patients from confined spaces or navigate stairways. The **TitanChair™** is designed to securely carry 375 lbs.




Taylor
TITAN
Series

TITANPC™



TITANPCX™



TITANCHAIR™



Item#	Product	Color	Material	Size	Packaging
43-TITANPC	TitanPC™ Soft Stretcher	White	Nonwoven/Poly	40" x 80"	10 per case
43-TITANPCX	TitanPCX™ Soft Stretcher	Safety Orange	Nonwoven/Poly	40" x 80"	10 per case
43-TITANCHAIR	TitanChair™	White	Nonwoven/Poly	36" x 36"	10 per case

*Fluid Impervious - Fluid does not penetrate the material even with time and pressure

**Backboard not included



LP43A0124