

SureFit™

DISPOSABLE FITTED STRETCHER SHEETS

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's durable SureFit™ fitted stretcher sheets are the best disposable fitted sheets available! Specially designed for patient transfer, Taylor fitted sheets provide the surest fit, gripping the mattress tightly and preventing slippage.

- Available in fluid resistant* and fluid impervious** versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets
- Our selection includes fitted sheets designed to aid in patient transfers up to 350 pounds
- Options available for standard and large cots, stretchers and gurneys

Taylor™

Healthcare Products, Inc.



Item #	Product	Color	Material	Size	Packaging
90-GFRC3690	G-Force™ Fitted Sheet, for Stryker® XPS	Blue	Fluid Resistant	36" x 90"	30/case
90-IFRC3690	I-Force™ Fitted Sheet, for Stryker® XPS	Blue	Fluid Impervious	36" x 90"	30/case
90-SUPER-G	SUPER-G™ Fitted Sheet, for Stryker® XPS	White	Fluid Resistant	36" x 90"	30/case
90-BSS3084	Standard Fitted Sheet 30"	Blue	Fluid Resistant	30" x 84"	50/case
90-BIS3084	Impervious Fitted Sheet 30"	Blue	Fluid Impervious	30" x 84"	50/case
90-BSS3484	Standard Fitted Sheet 34"	Blue	Fluid Resistant	34" x 84"	50/case
90-BIS3484	Impervious Fitted Sheet 34"	Blue	Fluid Impervious	34" x 84"	50/case
90-BXS3490	Extended Length Fitted Sheet	Blue	Fluid Resistant	34" x 90"	50/case
90-IXS3490	Impervious Extended Length Fitted Sheet	Blue	Fluid Impervious	34" x 90"	50/case
90-WXH3490	Heavy Duty Extended Length Fitted Sheet	White	Fluid Resistant	34" x 90"	50/case
90-WHS3484	Heavy Duty Fitted Sheet	White	Fluid Resistant	34" x 84"	50/case
90-IHS3484	Impervious Heavy Duty Fitted Sheet	White	Fluid Impervious	34" x 84"	30/case
92-SFS3084	Standard Fitted Sheet 250	Blue	Fluid Resistant	30" x 84"	50/case
92-IFS3084	Impervious Fitted Sheet 250	Blue	Fluid Impervious	30" x 84"	50/case
90-WFS4894	ER Fitted Sheet	White	Fluid Resistant	48" x 94"	50/case
90-BFS4894	ER Fitted Sheet	Blue	Fluid Resistant	48" x 94"	30/case

*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure

**Fluid Impervious - Fluid does not penetrate the material even with time and pressure

LP905FFSA0417