

SureFit™

**DISPOSABLE
FITTED
STRETCHER
SHEETS**



**Color Coded
Lift Capacity™**

Taylor®

Healthcare Products, Inc.

Taylor Healthcare Products, Inc. continues as the industry leader in quality, innovation and customer service!

Taylor Healthcare's durable SureFit™ fitted stretcher sheets are the best disposable fitted sheets available! Specially designed for patient transfer, Taylor fitted sheets provide the surest fit, gripping the mattress tightly and preventing slippage.

- The Color Coded sheets are designed to aid in patient transfers for the following weights when using proper lifting method:
 - Green – 250 lbs. Blue – 400 lbs.
 - Grey – 300 lbs. White – 450 lbs.
- Available in fluid resistant* and fluid impervious** versions
- Protects the mattress and patient against cross-contamination
- Elastic fitted design eliminates sheet slippage
- Soft linen-like feel maximizes patient comfort
- Reduces the overhead cost associated with reusable sheets
- Our selection includes fitted sheets designed to aid in patient transfers up to 450 pounds
- Options available for standard and large cots, stretchers & gurneys



Item #	Product	Color	Material	Size	Packaging
95-ECOFORCE	ECOFORCE™ Fitted Sheet for Stryker® XPS and Ferno® Power X1	Green	Fluid Resistant	36" x 90"	50/case
92-EFORCE	EFORCE™ Fitted Sheet for Stryker® XPS and Ferno® Power X1	Grey	Fluid Resistant	36" x 90"	50/case
90-GFRC3690	G-Force™ Fitted Sheet for Stryker® XPS and Ferno® Power X1	Blue	Fluid Resistant	36" x 90"	30/case
90-IFRC3690	I-Force™ Fitted Sheet for Stryker® XPS and Ferno® Power X1	Blue	Fluid Impervious	36" x 90"	30/case
90-SUPER-G	SUPER-G™ Fitted Sheet for Stryker® XPS and Ferno® Power X1	White	Fluid Resistant	36" x 90"	30/case
90-BSS3084	Standard Fitted Sheet 30"	Blue	Fluid Resistant	30" x 84"	50/case
90-BIS3084	Impervious Fitted Sheet 30"	Blue	Fluid Impervious	30" x 84"	50/case
90-BSS3484	Standard Fitted Sheet 34"	Blue	Fluid Resistant	34" x 84"	50/case
90-BIS3484	Impervious Fitted Sheet 34"	Blue	Fluid Impervious	34" x 84"	50/case
90-BXS3490	Extended Length Fitted Sheet	Blue	Fluid Resistant	34" x 90"	50/case
90-IXS3490	Impervious Extended Length Fitted Sheet	Blue	Fluid Impervious	34" x 90"	50/case
90-WXH3490	Heavy Duty Extended Length Fitted Sheet	White	Fluid Resistant	34" x 90"	50/case
90-WHS3484	Heavy Duty Fitted Sheet	White	Fluid Resistant	34" x 84"	50/case
90-IHS3484	Impervious Heavy Duty Fitted Sheet	White	Fluid Impervious	34" x 84"	30/case
92-SFS3084	Standard Fitted Sheet	Blue	Fluid Resistant	30" x 84"	50/case
92-IFS3084	Impervious Fitted Sheet	Blue	Fluid Impervious	30" x 84"	50/case
95-EFS3084	Economy Fitted Sheet	Green	Fluid Resistant	30" x 84"	50/case

*Fluid Resistant - Fluid beads up and could penetrate the material with time and pressure

**Fluid Impervious - Fluid does not penetrate the material even with time and pressure

LP90SFFSA0419