

# TITAN Body Bags

Taylor Healthcare Products, Inc. continues to be the industry leader in quality innovation and customer service!


Taylor Titan Body Bags set the new standard for efficiency, containment, functionality and strength. Our body bags have an identification pouch and consist of a six-handle design for durability during transportation.

- The Color Coded bags are designed and tested to aid in transportation for the following weights:

Black: 500 lbs. - 7 ft.

Royal Blue: 575 lbs. - 8 ft.

- Reduces risk of cross-contamination
- Fluid impervious\*
- Made of ultra-strong PVC for longevity
- An alternative-use as a cooling bag in heat-related situations



Item#	Product	Color	Material	Size	Packaging
46-BB500B	Titan™ Body Bag	Black	PVC	36" x 87"	5 per case
46-BB575RB	Titan™ Body Bag	Royal	Blue PVC	44" x 95"	5 per case

\* Fluid Impervious - Fluid does not penetrate the material even with time and pressure with the exception of the zipper