

Taylor Healthcare's Storage Trays attach to cots and provide a convenient, easy-to-reach area to place medical bags or other small, often used equipment. Storage trays keep your hands free to lift the stretcher, maintaining your back's health and improving patient care.

- Adds convenient storage space to your cot
- Constructed of durable nylon
- Secures easily to the cot with strong Velcro® fasteners
- Straps help secure equipment
- Head end trays come flat or with four storage pouches
- Fit most Stryker® cots



Item#	Product	Color	Material	Size	Packaging
45-BHT1010F	Flat Head End Tray	Black	Nylon	21" x 20"	10/case
45-BHT1010P	Head End Tray with Pouches	Black	Nylon	21" x 20"	10/case
45-BBT2010	Lower Storage Tray	Black	Nylon	38" x 18"	10/case

